



Molecular  
Dimensions

ACHIEVE MORE.

# Durham Osmolyte Screen

## Harness the power of intracellular environment-control compounds to improve your results.

The new Durham osmolyte screen allows you to test which of a broad range of cell components help stabilise your protein, improving the results of structural and functional studies.



Osmolytes are uncharged, highly soluble organic molecules that act in cells to balance high external salt concentrations without disturbing intracellular ionic strength. This makes osmolytes the ideal additives for stabilising your protein during structural and functional studies as they are usually present in high concentrations in native cell environments without affecting protein behaviour.

Easily identify which osmolytes are most beneficial to your research with The Durham Osmolyte screen, a 96-well screen that scans more than 50 different osmolytes, exclusively from Molecular Dimensions.

To optimise your protein's complete environment, we recommend using the Osmolyte screen together with the Durham pH and Salt screens. All three screens were developed at Durham University by Ehmke Pohl, Daniel Bruce and Emily Cardew.

### Clear presentation of results

Durham screen data can easily be analysed with NAMI. This free-to-download GUI-based software provides an intuitive plate-based colour gradient to enable easy identification of stabilising and denaturing conditions.

### Discover the detail

Bruce, D *et al.* *J. Vis. Exp* 144: e58666 (2019). Video presentation describing the use of stability screens in aiding your protein research. [Moldim.com/stability-video](https://moldim.com/stability-video).

### Osmolytes:

- > Often stabilise proteins as a result of their frequent presence in the native protein environment.
- > Rarely affect protein function because of their neutral chemical characteristics.
- > Uncharged, highly soluble organic molecules which microorganisms accumulate at significant concentrations to balance high salt external environments.
- > Include:
  - > Sugars
  - > Amino-acids
  - > Polyols
  - > NDSBs




Fig. 1: Examples of some of the osmolytes in the screen.

## ORDERING INFORMATION

	Pack Size	Description
MD1-120	96 x 0.5 mL	The Durham Osmolyte Screen
MD1-101	96 x 0.5 mL	The Durham pH Screen
MD1-102	96 X 0.5 mL	The Durham Salt Screen
MD1-103	384 x 0.5 mL	The Protein Stability Combo Kit (The Durham pH and Salt Screens and Rubik Buffer and Additive Screens)



Molecular  
Dimensions

ACHIEVE MORE.

# The Durham Salt and pH Screens

## Optimise your protein sample for better results from initial purification to final structure.

### How thermal shift assays help your research

By optimising protein stability you can improve your protein's crystallisability and the quality of your final structure. Thermal shift assays provide information on the effect various environmental factors have on the thermal denaturation of a protein, which in turn is related to protein stability. The most common thermal shift assay is ThermoFluor®, which measures increasing fluorescence induced by SYPRO® Orange binding to the exposed protein hydrophobic core.

### Streamline structure solution.

Salts are both common precipitating agents and additives in crystallisation conditions. Use the Durham Salt screen to identify which salts interact with your protein, and streamline the process of identifying the optimum crystallisation conditions.

### Ideal for fragment-based drug design.


The Durham pH Screen allows you to quickly identify molecular which behave like ligands and bind to specific pockets on the protein, competing with the small molecules in your fragment screen. Like ligands, these buffer molecules (de)stabilise your protein and can be detected by ThermoFluor.

	1	2	3	4	5	6	7	8	9	10	11	12
A	67.3	64.4	65.1	60.7	58.5	63.7	60.1	68.1	68.5	67.3	67.7	66.9
B	65.4	64.8	69.5	68.3	74.1	71.5	70.7	65.5	62.7	64.4	63.6	48.1
C	50.0	59.8	64.7	67.2	68.4	68.2	68.0	68.2	67.7	67.9	67.7	66.7
D	N/S	47.4	59.4	65.8	68.1	68.6	68.1	68.5	68.5	68.6	69.0	68.1
E	56.4	63.6	65.5	66.4	66.5	66.2	66.0	65.6	67.3	65.5	65.4	64.5
F	67.5	67.4	67.3	67.3	67.2	67.1	67.5	67.1	67.3	67.1	67.2	67.3
G	67.9	68.4	70.5	61.6	68.1	69.7	68.1	68.7	69.8	67.1	67.2	68.8
H	57.7	68.3	70.1	60.4	67.7	69.7	64.4	65.5	66.9	63.6	66.0	65.3

Least stable      No change      Most stable

Glucose Isomerase pH dependence

Dependence of melting temperature on concentration of magnesium and cobalt chloride for glucose isomerase. Thermal denaturation curves at a range of concentrations for (a)  $\text{CoCl}_2$  and (b)  $\text{MgCl}_2$ .



### Identify specifically bound metals and anions!

Using NAMI\* with the Salt Screen is a powerful tool for assessing the optimum concentration of salt to add to your protein. NAMI can generate waterfall plots of the thermal shift data at multiple salt concentrations or plots of melting temperature against concentration. These differentiate a sudden increase in protein stability indicating a specific interaction and steady increases in stability due to the changing protein environment.

\*NAMI is a free to download GUI-based program that aids thermal shift data analysis.



Part of Europa Biosite

### Your distributor in Switzerland

LubioScience GmbH  
Baumackerstrasse 24  
8050 Zurich

+41 41 417 02 80  
info@lubio.ch  
www.lubio.ch

ThermoFluor® is a registered trademark of Johnson & Johnson. SYPRO® Orange is a registered trademark of Life Technologies

A: Molecular Dimensions Inc  
434 West Dussel Drive  
Maumee OH 43537  
USA

T: +1 877 479 4339

F: +1 321 972-8896

W: moleculardimensions.com

E: enquiries@moleculardimensions.com