Silicon Carbide as a Novel RNA Affinity Medium with Improved Sensitivity and Size Diversity

Bernard Lam¹, Mohamed Elmogy², Song-Song Geng¹, Pam Roberts¹, Nezar Rghei¹ and Yousef Haj-Ahmad¹². ¹Norgen Biotek Corp., Thorold, ON, CANADA, ²Brock University, St. Catharines, ON, CANADA

Gene expression analysis using RNA has played a central role in our understanding of many aspects of the biology of an organism, including developmental timing, stem cell differentiation, signal transduction, disease and oncogenesis. In the past decade, studies of post-transcriptional gene regulation have brought attention to new classes of small RNAs, including microRNA (miRNA) and small interfering RNA (siRNA). The importance of small RNAs is clearly demonstrated by their diverse functions and applications. In fact, the emergence of the field of RNA interference has led to a drastic change in the criteria for RNA sample preparation methodologies. It is now highly important to be able to isolate total RNA with all sizes of RNA, including miRNA. Interestinor, many compressible variable trait RNA extraction products a charginal or and control or an extraction and control or an extraction or and the application of the control or an extraction or an extraction or control and control or an extraction or and the capital and the control of the co Gene expression analysis using RNA has played a central role in our understanding of Interestingly, many commercially available total RNA extraction products, particularly those using silica as a separation medium, selectively exclude RNA smaller than 200 those using silica as a separation medium, selectively exclude RNA smaller than 200 uncleatides. Additional steps, resulting in lengthy procedures, have to be incorporated to capture this important small RNA fraction. Here, we demonstrated that silicon carbide (SiC) can be utilized as a chromatographic medium for RNA purification. We showed that with the ues of SiC all sizes of RNA, including small RNA, can be isolated reproducibly with great sensitivity and ease. The true total RNA stoated using SiC, without the use of phase separation, consisted of similar quantity of miRNA to other miRNA isolation products utilizing a guanidium/phenol protocol. In addition, RNA (mRNA) and small miRNA when compared to a silica-based protocol or a traditional phase separation protocol. As SiC can be used in both a single column and high-throughput formats, it presents a robust solution to true total RNA isolation, particularly for rapid and simultaneous purification of miRNAs and their respective mRNA for rapid and simultaneous purification of miRNAs and their respective mRNA

В

600 -

Figure 2. Total RNA isolated using SiC Consists of a Similar Quantity and Diversity of Small RNA Species Compared to miRNA-Specific Products. Total RNA was isolated from 1 x 10° Heta cells using Norper's silicon cardide (Nor Total), or silica fiber (Supplier A Total) as the separation medium. Small RNA recovery was compared to two commercially available miRNA enrichment kits from the same suppliers (Nor miRNA and Supplier A miRNA). Guanidine thiocyanate/phenol-based TRI reagent isolation was used as a control for isolating the complete size range of RNA. Equal portions of isolated RNA were resolved on (A) a 1.5% formaldehyderagrose get or (B) an 9% Urena-PAGE gel. Only total RNA isolated using SiC (Panel A Lanes 1 and 2, Panel B Lanes 1 and 2) had the same amount of small RNA species as the kits that were specification general converts the nutrification and enrichment of



In the past decade, the emergence of the field of RNA interference has led to a drastic In the past decade, the emergence of the field of RNA interference has led to a drastic change in the criteria for RNA sample preparation methodologies. All sizes of RNA, including microRNA, should be able to be isolated with ease from a sample such that the expression level of both small RNAs and their mRNA targets can be correlated. The traditional method for total RNA isolation is the acid phenol/ guanidine thiocyanate method (Chomczynski and Sacchi, 1987), commonly known today as the Trizol (invitorgen, Carlsbad, CA) or TRI Reagent (Sigma-Aldrich, St. Louis, MO) method. This method delivers high yields of RNA and retrieves both small and large RNA species, however the method does rely on the use of harmful organic solvents and is rather time-consuming. During the last decade, there have been an increasing number of commerciality available siliculased since journm products for total RNA isolation. of commercially available silica-based spin column products for total RNA isolation. The isolation of RNA using these kits is rapid and generally does not require the use of organic solvents. However, the sizes of the RNA isolated with these kits are larger of organic solvents. However, the sizes of the RNA isolated with these kits are larger than 200 nuclotides, as most of the small RNA species, including miRNA and siRNA, are selectively excluded. Additional steps, such as incorporation of an acid phenol extraction step, are required for the capturing of small RNA molecules onto the silica media. Thus an increasing demand exists for a novel RNA isolation technology that can combine the species diversity of Trizol/TRI Reagent and the ease of use of spin column chromatography.

Objectives of the Study

- To demonstrate the use of silicon carbide (SiC) as an affinity medium for the purification of total RNA, including microRNAs, without use of phenol extraction.
 To compare the performance of SiC in RNA purification to that of current
- technology including phenol-based extraction and silica-based extraction.

Results






Figure 1. True Total RNA Isolation using SiC, without Phenol-Extraction. Total RNA was isolated from 1 x 10" HeLa cells using RNA isolation kits with spin column formats employing either silicon carbide (SiC), silica fiber (SiF) or silica matrix (SiM) as the separation medium. Guanidine thiocyanate/phenol-based TRI reagent (SiM) as the separation medium. Guanidine thicoyanate/phenol-based TRI reagent isolation was used as positive control for isolating the complete size range of RNA. Equal portions of isolated RNA were resolved on (A) a 1.5% formaldehyde-agarose glor or (B) an 8% Urea-PAGE gel. All kits isolated similar amounts of RNA as suggested by the intensity of the major ribosomal RNA bands. Only total RNA isolated using SiC (Panel A Lanes 1.2) or TRI reagent (TRI) (Panel A Lane 13, Panel B Lane 9) contained both the large and small RNA species, including those less than 200 n. In contrast to SiC, the RNA isolated using slice (Panel A Lanes 4-12, Panel B Lanes 3-8) included very few RNA species less than 200 nt in size, even though all kits compared employed similar lysis conditions based on Chomorzynski and Sacchi (1987). Nor = Norgen Total RNA Isolation Kit; TRI = Sigma-Aldrich TRI Reagent




Figure 3. High Recovery of Specific microRNA in Total RNA Isolated using Figure 3. High Recovery of Specific microRNA in Total RNA Isolated using SIC. The relative abundance of specific microRNA in the total RNA isolated using various kits was compared. RNA isolated in the aforementioned comparison (Figure 1) was polyaden/yated and used in RT-qPCR reactions using primers specific for the human miR-21 (Panel A), and miR-19 (Panel B). The resulting threshold cycle (C,) values were summarized in Panel D. Total RNA isolated using SiC contained levels of microRNA that are comparable to the levels obtained when RNA is isolated using TR. Reagent, as reflected by the C, values for the miR-21 and miR-19. In contrast, the C, values for the amplification of miR-21 and miR-19 using silica-bupfied total RNA are in the samplified miRNA of 515 at the same rate and gave similar C, values (Panel C), Δ = Silicon carbide-based Norgen Total RNA kit; Δ = Silica-based total RNA kits (Suppliers O, E and A); Δ = Sigma-Aldrich TRI Reagent. Curves in matching lighter colors represented corresponding –RT reactions and Δ = No Template Control





Contact Information:

Norgen Biotek 3430 Schmon Parkway Thorold, ON CANADA L2V 4Y6 Phone: (905) 227-8848 Fax: (905) 227-1061 nezar.rghei@norgenbiotek.com




Figure 4. Recovery of Specific microRNA in Total RNA Isolated using SiC was Comparable to that Isolated using microRNA Enrichment Products. The relative abundance of specific microRNA in the total RNA isolated with SiC was compared to that from miRNA enrichment products. RNA isolated in the aforementioned comparison (Figure 2) was polyadenlyated and used in RT-qPCR reactions using primers specific for the human miR-21 (Panel A), and miR-19 (Panel B). As indicated by the C, values (Panel D), both the microRNA purification kits and the SiC-based total RNA kit recovered microRNAs (both miR-19 and miR-21) that were very close to the levels obtained using TRI reagent. For the larger control SS rRNA, all the different separation methods showed amplification with similar C, values. Combining these results with the results in Figure 3, it can be seen that purification using SiC provides a new option for isolating a complete profile of RNA, including small microRNA. SiC allows for the simultaneous isolation and characterization of important microRNAs and their specific messenger RNA targets, without involving any phenol-based phase separation. ▲ "Silicon-based Norgen Total RNA Kit.

■ Norgen microRNA purification Kit.

■ Silican-based microRNA purification kit from Supplier A.

■ Sigma-Aldrich TRI
Reagent. Curves in matching lighter colors represented corresponding—RT reactions agent. Curves in matching lighter colors represented corresponding –RT reactions

| ____ = No Template Control




Figure 5. Total RNA Isolated using SiC Showed Better Linearity and Higher Sensitivity for both Large Messenger RNAs and Small microRNAs. Total RNA was isolated from increasing amounts of HeLa cells using a silicon carbide-based total RNA purification kit, a silica-fiber based total RNA purification kit, and guandine thiocyanatelyphenol-based TRI reagent as a control. The purified total RNA were then used in RT-Iq-CR to detect the microRNA miR-Z (Panel A) and the mRNA ST5 (Panel B). Silicon carbide columns were able to isolate both higher and more consistent levels of miR-Z transcripts over a wide input range of cells, down to as little as 10 cells (R² values: SiC = 0.9755 >> TRI = 0.9402 >> SiF = 0.4759). All three total RNA samples gave liner results for large messenger RNA species (human ST5 transcript tested), indicating the consistency of RNA isolation for all the three methods. Importantly, SiC-isolated total RNA again showed both the lowest C, values, as well as the highest R² value. \triangleq = Silicon carbide-based Norgen Total RNA Kit; \triangleq = Silicabased commercial total RNA kits (Supplier A); \triangleq = Sigma-Aldrich TRI Reagent.

Conclusions

- SiC-based RNA purification effectively isolated a complete profile of RNA species, from large messenger RNA to small microRNAs, without the use of phenol-extraction
- SiC-based RNA purification isolated both mRNA and miRNA with higher sensitivity as well as better linearity when compared to existing RNA purification technologies.

References

Chen, C. et al. 2005. Real-time quantification of microRNAs by stem-loop RT-PCR Nucleic Acids Research 33(20):e179

Chomczynski, P. & Sacchi, N. 1987. Single-step method of RNA isolation by acid guanidinium thiocyanate-phenol-chloroform extraction. Anal. Biochem. 162:156–159. Zhang, B., et al. 2007. microRNAs as oncogenes and tumor suppressors. Developmental Biology 302 (2007) 1–12.