

# efficiency in action



## Neon® Transfection System

The Neon® Transfection System is the next-generation electroporation system, designed to efficiently deliver DNA, RNA and protein into all cell types, especially difficult to transfect cell types, like stem cells and primary cells.



- Higher transfection efficiencies and viability in many cell types
- Greater flexibility— two tip sizes allow to transfect from  $2 \times 10^4$  up to  $1 \times 10^7$  cells per reaction  
Single universal reagent kit for all cell types, with 12-month shelf life
- Open system— tailor electroporation parameters to give the best results with your cells of interest
- Ease and speed of use thanks to the proprietary pipette chamber design

Product	Quantity	Cat. No.	List Price CHF	List Price CHF
Neon Transfection System	Each	MPK5000	6'370.00 CHF	5'414.50 CHF
Neon Transfection System Starter Pack	1 pack	MPK50000S	8'290.00 CHF	7'046.50 CHF



Life Technologies™ benchtop devices  
simple solutions to everyday complexities  
[www.lifetechnologies.com/benchtop](http://www.lifetechnologies.com/benchtop)

**15% off** list price\* on the Neon® Transfection System until 24/03/2012

Quote **P827643** when ordering