

BenchPro[®] 4100

Western Processing System

western detection
without the fuss

invitrogen[™]
by *life* technologies[™]



BenchPro[®] 4100 Western Processing System

The BenchPro[®] 4100 Western Processing System is designed to eliminate manual processing of western blot membranes.

By delivering consistent amounts of solution at precise times, the BenchPro[®] 4100 ensures reproducibility between experiments and eliminates possibility for errors. Finally, western processing made easy.

The system is programmable to execute any western protocol and each membrane can be processed with different set of reagents. The system is capable of processing from one to eight western membranes in parallel.



- Eliminates tedious hands-on work
- Provides more consistent results
- Eliminates possibility for errors
- Prevents contamination from one blot to the next
- Capable of executing any western detection protocol

Get **15% off** the BenchPro[®] 4100 until 24th March 2012 when you quote promo code **P827652** when ordering.



Life Technologies[™] benchtop devices
simple solutions to everyday complexities
www.lifetechnologies.com/benchtop

Product	Quantity	Cat. No.	List Price CHF	Promo Price CHF
BenchPro [®] 4100 Card Processing Station	Each	WP0001	7,590.00	6,451.50

lubio
science

To take advantage of this promotion, contact LuBioSciences info@lubio.ch