

15-1040-2/15-1070-2/15-1100-2

Cell strainers (40 µm/70 µm/100 µm diameter)

Size	100 pcs/box
Manufacturer	Biologix
Storage	Store at room temperature.
Shelf life	At least 12 months when stored properly.

Product Description

Cell strainers are sterile, easy to use devices for rapidly isolating primary cells to consistently obtain a uniform single-cell suspension from tissues. Fits 50 ml Falcon® conical tubes.

- Made of sturdy polypropylene and strong nylon mesh.
- Three different pore sizes available: 40µm, 70µm, or 100µm, strainers are color-coded according to pore size for convenience.
- Fits perfectly into most 50ml conical centrifuge tubes.
- Gamma radiation sterilized, non-pyrogenic.

Cat-No.	Size	Color	Pack size
15-1040-2	40µm	Purple	100 pcs/box
15-1070-2	70µm	Orange	100 pcs/box
15-1100-2	100µm	Green	100 pcs/box