

15-1040-2/15-1070-2/15-1100-2

Cell strainers (40 µm/70 µm/100 µm diameter)

| | |
|---------------------|--|
| Size | 100 pcs/box |
| Manufacturer | Biologix |
| Storage | Store at room temperature. |
| Shelf life | At least 12 months when stored properly. |

Product Description

Cell strainers are sterile, easy to use devices for rapidly isolating primary cells to consistently obtain a uniform single-cell suspension from tissues. Fits 50 ml Falcon® conical tubes.

- Made of sturdy polypropylene and strong nylon mesh.
- Three different pore sizes available: 40µm, 70µm, or 100µm, strainers are color-coded according to pore size for convenience.
- Fits perfectly into most 50ml conical centrifuge tubes.
- Gamma radiation sterilized, non-pyrogenic.

| Cat-No. | Size | Color | Pack size |
|-----------|-------|--------|-------------|
| 15-1040-2 | 40µm | Purple | 100 pcs/box |
| 15-1070-2 | 70µm | Orange | 100 pcs/box |
| 15-1100-2 | 100µm | Green | 100 pcs/box |