

Eligible Products Cusabio Promotion: Buy 2+ Cytokine ELISA Kits & Get 20% off

Code	Product Name	Species	Abbreviation	Sensitivity	Detection Range
CSB-EL001306HU	Human A disintegrin and metalloproteinase with thrombospondin motifs 18(ADAMTS18) ELISA kit	Human	ADAMTS18	0.016ug/ml	0.016ug/ml- 10ug/ml
CSB-E11848h	Human ADAM metallopeptidase with thrombospondin type 1 motif4,ADAMTS4 ELISA Kit	Human	ADAMTS4	3.125 ng/ml	12.5 ng/ml - 800 ng/ml
CSB-EL004774HU	Human C-C motif chemokine 1(CCL1) ELISA kit	Human	CCL1	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E14298m	Mouse chemokine (C-C motif) ligand 1 ELISA Kit	Mouse	CCL1	1.95 pg/ml	7.8 pg/ml - 500 pg/ml
CSB-E09762h	Human monocyte chemotactic protein 4,MCP-4 ELISA kit	Human	CCL13	9.75 pg/ml	39 pg/ml-2500 pg/ml
CSB-EL004777HU	Human C-C motif chemokine 14(CCL14) ELISA kit	Human	CCL14	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-EL004779HU	Human C-C motif chemokine 16(CCL16) ELISA kit	Human	CCL16	9.77 pg/mL	39.06 pg/mL-2500 pg/mL
CSB-E09257h	Human thymus activation regulated chemokine,TARC ELISA Kit	Human	CCL17	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E09941h	Human pulmonary activation regulated chemokine,PARC ELISA Kit	Human	CCL18	7.81 pg/ml	31.25 pg/ml-2000 pg/ml
CSB-E04671m	Mouse macrophage inflammatory protein 3β,MIP-3β ELISA kit	Mouse	CCL19	9.75 pg/ml	39 pg/ml - 2500 pg/ml
CSB-E04668m	Mouse Macrophage Inflammatory Protein 3α,MIP-3α ELISA kit	Mouse	CCL20	9.75 pg/ml	39 pg/ml - 2500 pg/ml
CSB-E04669r	Rat macrophage inflammatory protein 3α,MIP-3α ELISA Kit	Rat	CCL20	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E04660h	Human Macrophage-Derived Chemokine (MDC/CCL22) ELISA kit	Human	CCL22	10.324 pg/mL	15.625 pg/mL-1000 pg/mL
CSB-E04661m	Mouse Macrophage-Derived Chemokine,MDC ELISA kit	Mouse	CCL22	23.672 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E07278m	Mouse Eotaxin 2 ELISA Kit	Mouse	CCL24	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-EL004789MO	Mouse C-C motif chemokine 25(CCL25) ELISA kit	Mouse	CCL25	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-EL004790HU	Human C-C motif chemokine 26(CCL26) ELISA kit	Human	CCL26	9.75 pg/ml	39 pg/ml - 2500 pg/ml
CSB-E09125h	Human C-C motif chemokine 27 (CCL27) ELISA Kit	Human	CCL27	0.043 ng/mL	0.313 ng/mL-20 ng/mL
CSB-E07310h	Human mucosae associated epithelia chemokine,MEC ELISA Kit	Human	CCL28	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-EL004792MO	Mouse C-C motif chemokine 28(CCL28) ELISA kit	Mouse	CCL28	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E04662h	Human Macrophage Inflammatory Protein-1α,MIP-1α ELISA kit	Human	CCL3	9.75 pg/ml	39 pg/ml - 2500 pg/ml
CSB-E04663m	Mouse Macrophage Inflammatory Protein 1α,MIP-1α ELISA kit	Mouse	CCL3	1.95 pg/ml	7.8 pg/ml - 500 pg/ml
CSB-EL004797HU	Human C-C motif chemokine 4(CCL4) ELISA kit	Human	CCL4	7.8 pg/ml	31.2 pg/ml - 2000 pg/ml
CSB-E17375h	Human C-C motif chemokine 5 (CCL5/D17S136E/SCYA5) ELISA kit	Human	CCL5	2.307 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E17377h	Human C-C motif chemokine 7 (CCL7/MCP3/SCYA6/SCYA7) ELISA kit	Human	CCL7	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E04656h	Human monocyte chemotactic protein 2,MCP-2 ELISA kit	Human	CCL8	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-EL004936MO	Mouse Tumor necrosis factor receptor superfamily member 5(CD40) ELISA kit	Mouse	CD40	7.8 pg/ml	31.2 pg/ml - 2000 pg/ml
CSB-E05121h	Human decay-accelerating factor(DAF/CD55)ELISA Kit	Human	CD55	0.09 ng/ml	0.39 ng/ml - 25 ng/ml
CSB-E17547h	Human secretogranin-1 (chromogranin-B)(CHGB/SCG1)ELISA kit	Human	CHGB	0.39 ng/ml	1.56 ng/ml - 100 ng/ml
CSB-EL005381RA	Rat Muscarinic acetylcholine receptor M1(CHRM1) ELISA kit	Rat	CHRM1	0.27 ng/ml	0.27 ng/ml- 200ng/ml
CSB-E04659m	Mouse Macrophage Colony-Stimulating Factor,M-CSF ELISA kit	Mouse	CSF1	1.95 pg/mL	7.81 pg/mL-500 pg/mL
CSB-E04658h	Human Macrophage Colony-Stimulating Factor,M-CSF ELISA kit	Human	CSF1	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E07424r	Rat Macrophage Colony-Stimulating Factor,M-CSF ELISA kit	Rat	CSF1	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E10012h	Human Macrophage Colony-Stimulating Factor Receptor,M-CSFR ELISA Kit	Human	CSF1R	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E04568h	Human Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit	Human	CSF2	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E04569m	Mouse Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit	Mouse	CSF2	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E04570r	Rat Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit	Rat	CSF2	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E04564m	Mouse Granulocyte Colony Stimulating Factor,G-CSF ELISA Kit	Mouse	CSF3	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E04563h	Human Granulocyte Colony Stimulating Factor,G-CSF ELISA Kit	Human	CSF3	1.95 pg/ml	7.8 pg/ml - 500 pg/ml
CSB-E17503h	Human Cystatin-A(CSTA/STF1/STFA)ELISA kit	Human	CSTA	0.156 ng/ml	0.625 ng/ml- 40ng/ml
CSB-E09150h	Human growth-regulated oncogene/melanoma growth stimulating activity,GROα/MGSA ELISA Kit	Human	CXCL1	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E17286m	Mouse chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha) (CXCL1) ELISA kit	Mouse	CXCL1	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E08181h	Human interferon-inducible protein 10,IP-10 ELISA Kit	Human	CXCL10	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E08183m	Mouse interferon-inducible protein 10,IP-10 ELISA Kit	Mouse	CXCL10	0.078 ng/ml	0.312 ng/ml - 20 ng/ml
CSB-E09023h	Human Interferon inducible T-cell Chemoattractant,I-TAC ELISA Kit	Human	CXCL11	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E16832m	Mouse C-X-C motif chemokine 13 (Cxcl13/Blc/Scyb13) ELISA kit	Mouse	CXCL13	0.096 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E10025h	Human Breast and kidney expressed chemokine,BRAK ELISA Kit	Human	CXCL14	0.09 ng/ml	0.39 ng/ml-25 ng/ml
CSB-E08871h	Human CXC-chemokine ligand 16,CXCL16 ELISA Kit	Human	CXCL16	3.603 pg/mL	46.875 pg/mL-3000 pg/mL
CSB-E08873m	Mouse CXC-chemokine ligand 16,CXCL16 ELISA Kit	Mouse	CXCL16	3.012 pg/mL	7.8125 pg/mL-500 pg/mL
CSB-E04666m	Mouse Macrophage Inflammatory Protein 2,MIP-2 ELISA kit	Mouse	CXCL2	5.8 pg/mL	23.5 pg/mL-1500 pg/mL
CSB-EL006249RA	Rat C-X-C motif chemokine 3(CXCL3) ELISA kit	Rat	CXCL3	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E09024h	Human monocyte interferon gamma inducing factor,MIG ELISA Kit	Human	CXCL9	9.75 pg/ml	39 pg/ml-2500 pg/ml
CSB-EL006252MO	Mouse C-X-C motif chemokine 9(CXCL9) ELISA kit	Mouse	CXCL9	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E12825h	Human CXC-chemokine receptor 4(CXCR4)ELISA Kit	Human	CXCR4	0.39 pg/mL	1.56 pg/mL-100 pg/mL
CSB-E08027h	Human Epidermal growth factor,EGF ELISA Kit	Human	EGF	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E08028m	Mouse Epidermal growth factor,EGF ELISA Kit	Mouse	EGF	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E12124h	Human epidermal growth factor receptor,EGFR ELISA Kit	Human	EGFR	0.019 ng/mL	0.312 ng/mL-20 ng/mL
CSB-E04538h	Human Erythropoietin,EPO ELISA Kit	Human	EPO	Request Information	Request Information
CSB-E04537h	Human Erythropoietin receptor,EPOR ELISA Kit	Human	EPOR	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E04546h	Human acidic fibroblast growth factor,aFGF/FGF-1 ELISA Kit	Human	FGF1	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E07327m	Mouse acidic fibroblast growth factor,aFGF/FGF-1 ELISA Kit	Mouse	FGF1	1.56 pg/mL	6.25 pg/mL-400 pg/mL
CSB-EL522052MO	Mouse Fibroblast growth factor 15(FGF15) ELISA kit	Mouse	FGF15	3.9 pg/mL	15.6 pg/mL - 1000 pg/ml
CSB-EL008624HU	Human Fibroblast growth factor 19(FGF19) ELISA kit	Human	FGF19	6.25 pg/mL	25 pg/mL-1600 pg/mL
CSB-E08002r	Rat basic fibroblast growth factor,bFGF ELISA Kit	Rat	FGF2	0.078 ng/ml	0.312 ng/ml - 20 ng/ml
CSB-EL008627MO	Mouse Fibroblast growth factor 21(FGF21) ELISA kit	Mouse	FGF21	7.81 pg/ml	31.25 pg/ml - 2000 pg/ml
CSB-E12170r(1)	Rat fibroblast growth factor 23 (FGF23) ELISA Kit	Rat	FGF23	3.9 ng/mL	15.6 ng/mL-1000 ng/mL
CSB-E11859h	Human Vitamin D-binding protein,DBP ELISA Kit	Human	GC	0.078 ng/ml	0.3125 ng/ml - 20ng/ml
CSB-E12009h	Human growth differentiation factor 15(GDF15)ELISA Kit	Human	GDF15	1.95 pg/mL	7.8 pg/mL-500 pg/mL

CSB-EL009345MO	Mouse Growth/differentiation factor 15(GDF15) ELISA kit	Mouse	GDF15	1.95 pg/ml	7.8 pg/ml - 500 pg/ml
CSB-E11769h	Human growth differentiation factor 2,GDF2 ELISA Kit	Human	GDF2	0.98 pg/mL	3.9 pg/mL-250 pg/mL
CSB-EL009703HU	Human glypican 1 (GPC1) ELISA kit	Human	GPC1	0.039 ng/ml	0.156 ng/ml-10 ng/ml
CSB-E08160h	Human Gelsolin ELISA Kit	Human	GSN	2.74 ng/ml	2.74 ng/ml - 2000 ng/ml
CSB-E04573h	Human hepatocyte growth factor,HGF ELISA kit	Human	HGF	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E07347m	Mouse hepatocyte growth factor,HGF ELISA kit	Mouse	HGF	Request Information	Request Information
CSB-E07346r	Rat hepatocyte growth factor,HGF ELISA kit	Rat	HGF	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-EL010767HU	Human 3-hydroxyacyl-CoA dehydrogenase type-2(HSD17B10) ELISA kit	Human	HSD17B10	3.9 pg/ml	15.625 pg/ml- 1000pg/ml
CSB-EL011038HU	Human interferon alpha-2(IFNA2) ELISA kit	Human	IFNA2	0.78 pg/mL	3.12 pg/mL-200 pg/mL
CSB-EL011047HU	Human Interferon alpha/beta receptor 2(IFNAR2) ELISA kit	Human	IFNAR2	7.8 pg/ml	31.25 pg/ml - 2000pg/ml
CSB-E04578m(1)	Mouse Interferon γ ,IFN- γ ELISA Kit	Mouse	IFNG	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E04577h	Human Interferon γ (IFN- γ) ELISA Kit	Human	IFNG	0.078 ng/ml	0.312 ng/ml-20 ng/ml
CSB-E14210h	Human interferon γ (IFN- γ) antibody ELISA kit	Human	IFN- γ Ab	Request Information	Request Information
CSB-E04581m	Mouse Insulin-like growth factor 1,IGF-1 ELISA Kit	Mouse	IGF1	0.151 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E04582r	Rat Insulin-like growth factor 1,IGF-1 ELISA Kit	Rat	IGF1	0.151 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E04586h	Human insulin-like growth factors binding protein 1,IGFBP-1 ELISA Kit	Human	IGFBP1	0.39 ng/mL	1.56 ng/mL-100 ng/mL
CSB-E04588h	Human insulin-like growth factors binding protein 2,IGFBP-2 ELISA Kit	Human	IGFBP2	0.195 ng/mL	0.78 ng/mL-50 ng/mL
CSB-E04589m	Mouse Insulin-like growth factor binding protein 2,IGFBP-2 ELISA KIT	Mouse	IGFBP2	0.11 ng/mL	0.45 ng/mL-30 ng/mL
CSB-E07351r	Rat Insulin-like growth factor binding protein 2,IGFBP-2 ELISA kit	Rat	IGFBP2	0.156 ng/mL	0.625 ng/mL-400 ng/mL
CSB-E04590h	Human insulin-like growth factors binding protein 3,IGFBP-3 ELISA Kit	Human	IGFBP3	0.195 ng/mL	0.78 ng/mL-50 ng/mL
CSB-E04591m	Mouse Insulin-like growth factor binding protein 3,IGFBP-3 ELISA KIT	Mouse	IGFBP3	0.196 ng/mL	0.391 ng/mL-25 ng/mL
CSB-E07353r	Rat Insulin-like growth factor binding protein 3,IGFBP-3 ELISA kit	Rat	IGFBP3	2.7 ng/ml	2.7 ng/ml - 2000 ng/ml
CSB-E04592h	Human insulin-like growth factors binding protein 4,IGFBP-4 ELISA Kit	Human	IGFBP4	15.6 ng/mL	62.5 ng/mL-4000 ng/mL
CSB-EL011099MO	Mouse Insulin-like growth factor-binding protein 5(IGFBP5) ELISA kit	Mouse	IGFBP5	19.5 pg/ml	78 pg/ml - 5000pg/ml
CSB-E12784r	Rat Insulin-like growth factor binding protein 5,IGFBP-5 ELISA kit	Rat	IGFBP5	0.078 ng/mL	0.312 ng/mL-20 ng/mL
CSB-EL011100HU	Human Insulin-like growth factor-binding protein 6(IGFBP6) ELISA kit	Human	IGFBP6	0.195 ng/mL	0.78 ng/mL-50 ng/mL
CSB-E17249h	Human Insulin-like growth factor-binding protein 7 (IGFBP7/MAC25/PSF) ELISA kit	Human	IGFBP7	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E12015B	Bovine Immunoglobulin G,IgG ELISA Kit	Bovine	IgG	19.5 ng/ml	78 ng/ml-5000ng/ml
CSB-E04593h	Human Interleukin 10,IL-10 ELISA KIT	Human	IL10	7.8 pg/ml	31.25 pg/ml - 2000 pg/ml
CSB-E04594m	Mouse interleukin 10,IL-10 ELISA KIT	Mouse	IL10	15.6 pg/ml	62.5 pg/ml - 4000 pg/ml
CSB-E04597m	Mouse Interleukin 11,IL-11 ELISA KIT	Mouse	IL11	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E04600m	Mouse Interleukin 12,IL-12/P70 ELISA KIT	Mouse	IL12/P70	15.6 pg/ml	62.5 pg/ml - 4000 pg/ml
CSB-E13288p	Pig Interleukin 12,IL-12/P70 ELISA Kit	Pig	IL12/P70	0.11 ng/mL	0.45 ng/mL-30 ng/mL
CSB-E04599h	Human Interleukin 12,IL-12/P70 ELISA KIT	Human	IL-12/P70	7.8 pg/ml	31.25 pg/ml - 2000pg/ml
CSB-E04598h	Human Interleukin 12,IL-12/P40 ELISA KIT	Human	IL12B	0.039 ng/mL	0.156 ng/ml-10 ng/ml
CSB-E07360m	Mouse interleukin 12,IL-12/P40 ELISA KIT	Mouse	IL12B	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E04602m	Mouse Interleukin 13,IL-13 ELISA Kit	Mouse	IL13	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E04601h	Human Interleukin 13,IL-13 ELISA KIT	Human	IL13	0.09 ng/ml	0.39 ng/ml-25 ng/ml
CSB-E04603h	Human Interleukin 15,IL-15 ELISA KIT	Human	IL15	0.78 pg/ml	3.12 pg/ml - 200 pg/ml
CSB-E04604m	Mouse Interleukin 15,IL-15 ELISA Kit	Mouse	IL15	0.078 ng/ml	0.3125ng/ml - 20ng/ml
CSB-E04605h	Human Interleukin 16,IL-16 ELISA KIT	Human	IL16	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E04608m	Mouse Interleukin 17,IL-17 ELISA Kit	Mouse	IL17A	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E12819h	Human Interleukin 17A (IL-17A/IL-17) ELISA Kit	Human	IL17A	15.6 pg/ml	62.5 pg/ml - 4000 pg/ml
CSB-E15918h	Human Interleukin17B (IL17B) ELISA Kit	Human	IL17B	7.8 pg/mL	31.25 pg/ml-2000 pg/ml
CSB-E15917h	Human Interleukin17C (IL17C) ELISA Kit	Human	IL17C	0.039 ng/ml	0.156 ng/ml -10 ng/ml
CSB-EL011601MO	Mouse Interleukin-17F(IL17F) ELISA kit	Mouse	IL17F	7.81 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-EL011602HU	Human Interleukin-17 receptor A(IL17RA) ELISA kit	Human	IL17RA	9.75 pg/mL	39 pg/ml-2500 pg/ml
CSB-EL011603HU	Human Interleukin-17 receptor B(IL17RB) ELISA kit	Human	IL17RB	0.078 ng/ml	0.3125ng/ml- 20ng/ml
CSB-E07450h	Human Interleukin 18,IL-18 ELISA KIT	Human	IL18	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E13687h	Human interleukin-18 binding proreïn(IL-18BP) ELISA Kit	Human	IL18BP	3.9 pg/ml	15.6 pg/ml-1000 pg/ml
CSB-E04621m	Mouse Interleukin 1 α ,IL-1 α ELISA kit	Mouse	IL1A	0.58 pg/ml	2.35 pg/ml - 150 pg/ml
CSB-E04622r	Rat Interleukin 1 α ,IL-1 α ELISA kit	Rat	IL1A	19.5 pg/ml	78 pg/ml-5000 pg/ml
CSB-E08053h	Human Interleukin 1 β ,IL-1 β ELISA Kit	Human	IL1B	1.95 pg/ml	7.8 pg/ml - 500 pg/ml
CSB-E08055r	Rat Interleukin 1 β ,IL-1 β ELISA Kit	Rat	IL1B	0.039 ng/ml	0.156 ng/ml-10 ng/ml
CSB-EL011616HU	Human Interleukin-36 receptor antagonist protein(IL1F5) ELISA kit	Human	IL1F5	7.8 pg/ml	31.25 pg/ml - 2000 pg/ml
CSB-E04611h	Human soluble interleukin-1 receptor I ,IL-1sR I ELISA kit	Human	IL1R1	0.039 ng/ml	0.156 ng/ml-10 ng/ml
CSB-E04612h	Human soluble interleukin-1 receptor II ,IL-1sR II ELISA kit	Human	IL1R2	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-EL011626MO	Mouse Interleukin-1 receptor-like 1(IL1RL1) ELISA kit	Mouse	IL1RL1	19.5 pg/ml	78 pg/ml-5000 pg/ml
CSB-E10395m	Mouse Interleukin 1 receptor antagonist, IL-1ra ELISA kit	Mouse	IL1RN	15.6 pg/mL	62.5 pg/mL-4000 pg/mL
CSB-E10396h	Human Interleukin 1 receptor antagonist, IL-1ra ELISA kit	Human	IL1RN	3.9 pg/ml	15.6 pg/ml-1000 pg/ml
CSB-E04626h	Human Interleukin 2,IL-2 ELISA kit	Human	IL2	7.81 pg/ml	31.25 pg/ml-2000 pg/ml
CSB-E04627m	Mouse Interleukin 2,IL-2 ELISA kit	Mouse	IL2	7.8 pg/ml	31.25 pg/ml - 2000 pg/ml
CSB-E15015h	Human Interleukin 20(IL-20)ELISA Kit	Human	IL20	0.078 ng/ml	0.3125 ng/ml- 20ng/ml
CSB-EL011630MO	Mouse Interleukin-20(IL20) ELISA kit	Mouse	IL20	7.81 pg/ml	31.25 pg/ml-2000 pg/ml
CSB-E11707h	Human Interleukin 21,IL-21 ELISA Kit	Human	IL21	0.78 pg/mL	3.12 pg/mL-200 pg/mL
CSB-E13162m	Mouse Interleukin 21(IL-21) ELISA Kit	Mouse	IL21	15.6 pg/ml	62.5 pg/ml-4000 pg/ml
CSB-E13513m	Mouse Interleukin 22(IL-22) ELISA Kit	Mouse	IL22	0.078 ng/mL	0.312 ng/mL-20 ng/mL
CSB-E13418h	Human Interleukin 22(IL-22) ELISA KIT	Human	IL22	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E08461h	Human Interleukin 23,IL-23 ELISA Kit	Human	IL-23	1.56 pg/mL	6.25 pg/mL-400 pg/mL
CSB-E15840h	Human Interleukin 24(IL-24)ELISA Kit	Human	IL24	0.078 ng/ml	0.3125 ng/ml- 20ng/ml
CSB-E08464h	Human Interleukin 27,IL-27 ELISA Kit	Human	IL27	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E11708h	Human Interleukin 28A,IL-28A/IFN- λ 2 ELISA kit	Human	IL-28A/IFN- λ 2	0.039 ng/ml	0.156 ng/ml - 10 ng/ml

CSB-E14290h	Human Interleukin 29(IL-29)ELISA Kit	Human	IL29	15.6 pg/ml	62.5 pg/ml-4000 pg/ml
CSB-E04629h	Human soluble interleukin-2 receptor,IL-2sR α ELISA kit	Human	IL2RA	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-E07440m	Mouse soluble interleukin-2 receptor,sIL-2R α ELISA kit	Mouse	IL2RA	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E04631h	Human Interleukin 3,IL-3 ELISA KIT	Human	IL3	15.6 pg/ml	62.5 pg/ml-4000 pg/ml
CSB-E04632m	Mouse Interleukin 3,IL-3 ELISA KIT	Mouse	IL3	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E16999h	Human Interleukin 31(IL-31) ELISA Kit	Human	IL31	31.25 pg/ml	125 pg/ml-8000 pg/ml
CSB-E13000h	Human Interleukin 33(IL-33) ELISA Kit	Human	IL33	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E14393m	Mouse Interleukin 33(IL-33) ELISA Kit	Mouse	IL33	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E04633h	Human Interleukin 4,IL-4 ELISA KIT	Human	IL4	7.8 pg/ml	31.25 pg/ml - 2000 pg/ml
CSB-E04634m	Mouse Interleukin 4,IL-4 ELISA KIT	Mouse	IL4	7.8 pg/ml	31.25 pg/ml-2000 pg/ml
CSB-E04635r	Rat Interleukin 4,IL-4 ELISA KIT	Rat	IL4	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E04636h	Human Interleukin 5,IL-5 ELISA KIT	Human	IL5	0.97 pg/ml	3.9 pg/ml-250 pg/ml
CSB-E04637m	Mouse Interleukin 5,IL-5 ELISA KIT	Mouse	IL5	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E04638h	Human Interleukin 6,IL-6 ELISA KIT	Human	IL6	2.453 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E16524h	Human Interleukin 6 receptor(IL-6R)ELISA Kit	Human	IL6R	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E14032h	Human Interleukin-7,IL-7 ELISA kit	Human	IL7	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E10257m	Mouse Interleukin-7,IL-7 ELISA kit	Mouse	IL7	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-EL011670MO	Mouse Interleukin-7 receptor subunit alpha(IL7R) ELISA kit	Mouse	IL7R	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E08052c	Dog Interleukin-8,IL-8 ELISA Kit	Dog	IL8	0.146 ng/mL	0.312 ng/mL-20 ng/mL
CSB-E04641h	Human Interleukin 8,IL-8 ELISA KIT	Human	IL8	7.110 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E06787p	Pig interleukin 8,IL-8 ELISA Kit	Pig	IL8	31.25 pg/mL	125 pg/mL-8000 pg/mL
CSB-EQ027463HU	Human interleukin 8(IL8) autoantibody ELISA kit	Human	IL8 auto-Ab	Request Information	Request Information
CSB-EL012824HU	Human Plastin-2(LCP1) ELISA kit	Human	LCP1	3.9 pg/ml	15.6 pg/ml-1000 pg/ml
CSB-EL013042HU	Human Lysyl oxidase homolog 3(LOXL3) ELISA kit	Human	LOXL3	0.078 ng/ml	0.312ng/ml - 20ng/ml
CSB-EL013323HU	Human Melanoma-associated antigen 1(MAGEA1) ELISA kit	Human	MAGEA1	0.039 ng/ml	0.156 ng/ml-10ng/ml
CSB-E08283h	Human myelin basic protein,MBP ELISA Kit	Human	MBP	1.95 ng/ml	7.8 ng/ml- 500ng/ml
CSB-EL013572HU	Human Methylcrotonoyl-CoA carboxylase subunit alpha, mitochondrial(MCC1) ELISA kit	Human	MCC1	0.19 ng/ml	0.78 ng/ml-50ng/ml
CSB-E13490h	Human hepatocyte growth factor receptor(c-MET/HGFR)ELISA Kit	Human	MET	7.8 pg/mL	31.2 pg/mL-2000 pg/mL
CSB-E13492m	Mouse Hepatocyte Growth Factor Receptor(C-MET/HGFR)ELISA Kit	Mouse	MET	0.02 ng/mL	0.078 ng/mL-5 ng/mL
CSB-EL013750HU	Human Microfibrillar-associated protein 5(MFAP5) ELISA kit	Human	MFAP5	3.125pg/ml	12.5pg/ml - 800pg/ml
CSB-E04683h	Human Nerve growth factor,NGF ELISA kit	Human	NGF	6.86 pg/ml	6.86 pg/ml-5000 pg/ml
CSB-E04684m	Mouse Nerve growth factor,NGF ELISA kit	Mouse	NGF	15.6 pg/ml	62.5 pg/ml-4000 pg/ml
CSB-E09886h	Human Nephritin ELISA Kit	Human	NPHS1	0.39 ng/ml	1.56 ng/ml-100ng/ml
CSB-E09144h	Human Neurotensin,NT ELISA Kit	Human	NTS	3.9 pg/ml	15.6 pg/ml - 1000pg/ml
CSB-E09891h	Human Podocalyxin,PCX ELISA Kit	Human	PCX	0.54 ng/ml	0.54 ng/ml-400 ng/ml
CSB-EL017667HU	Human Programmed cell death 1 ligand 2(PDCD1LG2) ELISA kit	Human	PDCD1LG2	0.11 ng/ml	0.45 ng/ml-30 ng/ml
CSB-E14313r	Rat Platelet-Derived Growth Factor A,PDGF-A ELISA kit	Rat	PDGFA	7.8 pg/ml	31.25pg/ml-2000pg/ml
CSB-E04701h	Human Platelet-Derived Growth Factor AB,PDGF-AB ELISA kit	Human	PDGF-AB	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E08923h	Human Platelet-Derived Growth Factor-BB,PDGF-BB ELISA kit	Human	PDGF-BB	7.81 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E16491h	Human Profilin 1 ELISA Kit	Human	PFN1	0.195 ng/ml	0.78 ng/ml - 50 ng/ml
CSB-E15980h	Human Prolactin Inducible Protein(PIP)ELISA Kit	Human	PIP	0.055 ng/ml	0.055 ng/ml - 40ng/ml
CSB-E13155h	Human phospholipase A2 receptor 1, 180kDa (PLA2R1) ELISA kit	Human	PLA2R1	0.027 ng/ml	0.027 ng/ml - 20ng/ml
CSB-E17578h	Human eosinophil granule major basic protein (EMBP/MBP) ELISA Kit	Human	PRG2	1.3 ng/ml	1.3 ng/ml - 1000 ng/ml
CSB-EL018915HU	Human Proteasome activator complex subunit 1(PSME1) ELISA kit	Human	PSME1	0.049 ng/ml	0.19 ng/ml-12.5ng/ml
CSB-EL020808HU	Human Secretogranin-2(SCG2) ELISA kit	Human	SCG2	0.59 ng/ml	2.34 ng/ml - 150ng/ml
CSB-EL021135HU	Human 14-3-3 protein sigma(SFN) ELISA kit	Human	SFN	0.27 ng/ml	0.27 ng/ml- 200ng/ml
CSB-EL022411HU	Human Sortilin-related receptor(SORL1) ELISA kit	Human	SORL1	1.95 ng/ml	7.8 ng/ml- 500ng/ml
CSB-EL022575HU	Human Pancreatic secretory trypsin inhibitor(SPINK1) ELISA kit	Human	SPINK1	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E14144m	Mouse thymus activation regulated chemokine,TARC ELISA Kit	Mouse	TARC/CCL17	10.512 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-EL023131HU	Human Threonyl-tRNA synthetase, cytoplasmic(TARS) ELISA kit	Human	TARS	0.039 ng/ml	0.156ng/ml - 10 ng/ml
CSB-E04724h	Human transforming growth factor α ,TGF- α ELISA Kit	Human	TGFA	15.6 pg/mL	62.5 pg/mL-4000 pg/mL
CSB-E07289r	Rat transforming growth factor α ,TGF- α ELISA Kit	Rat	TGFA	1.95 pg/ml	7.8 pg/ml-500pg/ml
CSB-E04796c	Dog Transforming growth factor β 1,TGF- β 1 ELISA Kit	Dog	TGFB1	0.747 ng/ml	0.781 ng/ml - 50 ng/ml
CSB-E04725h	Human Transforming Growth factor β 1,TGF- β 1 ELISA kit	Human	TGFB1	0.747 ng/mL	0.78 ng/mL-50 ng/mL
CSB-E14208B	Bovine transforming growth factor β 1,TGF- β 1 ELISA Kit	Bovine	TGFB1	1.56 ng/mL	6.25 ng/mL-400 ng/mL
CSB-E04726m	Mouse Transforming Growth factor β 1,TGF- β 1 ELISA kit	Mouse	TGFB1	0.2 ng/mL	0.78 ng/mL-50 ng/mL
CSB-EQ027985HU	Human Transforming growth factor beta-1(TGFB1) autoantibody ELISA kit	Human	TGFB1 auto-Ab	Request Information	Request Information
CSB-E09783h	Human transforming growth factors β 2,TGFB2 ELISA Kit	Human	TGFB2	Request Information	Request Information
CSB-EQ027984HU	Human transforming growth factor, beta 2 (TGFB2) autoantibody ELISA kit	Human	TGFB2 auto-Ab	Request Information	Request Information
CSB-E16665h	Human Transforming growth factor-beta-induced protein ig-h3 (TGFB1/BIGH3) ELISA kit	Human	TGFB1	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-EL023450MO	Mouse Transforming growth factor-beta-induced protein ig-h3(TGFB1) ELISA kit	Mouse	TGFB1	Request Information	Request Information
CSB-E13669h	Human transforming growth factor β receptor 2 (TGF- β R2) ELISA kit	Human	TGFB2	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-EL023453HU	Human Transforming growth factor beta receptor type 3(TGFB3) ELISA kit	Human	TGFB3	Request Information	Request Information
CSB-E16980p	Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit	Pig	TNF	0.040 ng/mL	0.047 ng/mL-3 ng/mL
CSB-E04741m(1)	Mouse Tumor necrosis factor α ,TNF- α ELISA KIT	Mouse	TNF	15.6 pg/mL	62.5 pg/mL-4000 pg/mL
CSB-E04740h	Human Tumor necrosis factor α ,TNF- α ELISA KIT	Human	TNF	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-EL023959HU	Human Tumor necrosis factor-inducible gene 6 protein(TNFAIP6) ELISA kit	Human	TNFAIP6	11.75 ng/ml	47 ng/ml - 3000 ng/ml
CSB-E04749h	Human tumor necrosis factor-related apoptosis-inducing ligand receptor 4(TRAIL-R4) ELISA Kit	Human	TNFRSF10D	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E13539h	Human Soluble Receptor Activator of Nuclear Factor-kB,sRANK ELISA Kit	Human	TNFRSF11A	0.078 ng/ml	0.3125 ng/ml- 20ng/ml
CSB-EL023971HU	Human Tumor necrosis factor receptor superfamily member 13B(TNFRSF13B) ELISA kit	Human	TNFRSF13B	19.5 pg/ml	78 pg/ml - 5000 pg/ml

CSB-EL023971MO	Mouse Tumor necrosis factor receptor superfamily member 13B(TNFRSF13B) ELISA kit	Mouse	TNFRSF13B	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-EL023972MO	Mouse Tumor necrosis factor receptor superfamily member 13C(TNFRSF13C) ELISA kit	Mouse	TNFRSF13C	9.75 pg/ml	39 pg/ml - 2500 pg/ml
CSB-EL023973HU	Human Tumor necrosis factor receptor superfamily member 14(TNFRSF14) ELISA kit	Human	TNFRSF14	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-EL023974HU	Human Tumor necrosis factor receptor superfamily member 17(TNFRSF17) ELISA kit	Human	TNFRSF17	1.95 pg/ml	7.8 pg/ml- 500 pg/ml
CSB-EL023974MO	Mouse Tumor necrosis factor receptor superfamily member 17(TNFRSF17) ELISA kit	Mouse	TNFRSF17	0.09 ng/ml	0.39 ng/ml- 25ng/ml
CSB-EL023975HU	Human Tumor necrosis factor receptor superfamily member 18(TNFRSF18) ELISA kit	Human	TNFRSF18	0.039 ng/ml	0.156 ng/ml- 10 ng/ml
CSB-EL023975MO	Mouse Tumor necrosis factor receptor superfamily member 18(TNFRSF18) ELISA kit	Mouse	TNFRSF18	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E04736h	Human Tumor necrosis factor soluble receptor I ,TNFsR- I ELISA KIT	Human	TNFRSF1A	65 pg/mL	78.125 pg/mL-5000 pg/mL
CSB-E04737m	Mouse Tumor necrosis factor soluble receptor I ,TNFsR- I ELISA KIT	Mouse	TNFRSF1A	1.261 pg/mL	15.625 pg/mL-1000 pg/mL
CSB-E11266h	Human soluble tumor necrosis factor receptor 2,sTNF-R2 ELISA Kit	Human	TNFRSF1B	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E04739m	Mouse Tumor necrosis factor soluble receptor II ,TNFsR- II ELISA KIT	Mouse	TNFRSF1B	1.6 pg/mL	7.8125 pg/mL-500 pg/mL
CSB-EL023979HU	Human Tumor necrosis factor receptor superfamily member 21(TNFRSF21) ELISA kit	Human	TNFRSF21	9.75 pg/mL	39 pg/mL-2500 pg/mL
CSB-EL023982HU	Human Tumor necrosis factor receptor superfamily member 6B(TNFRSF6B) ELISA kit	Human	TNFRSF6B	0.156 ng/ml	0.625 ng/ml - 40 ng/ml
CSB-EL023984HU	Human Tumor necrosis factor receptor superfamily member 9(TNFRSF9) ELISA kit	Human	TNFRSF9	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-EL023984MO	Mouse Tumor necrosis factor receptor superfamily member 9(TNFRSF9) ELISA kit	Mouse	TNFRSF9	3.9 pg/ml	15.6 pg/ml - 1000 pg/ml
CSB-E04750h	Human soluble tumor necrosis factor-related apoptosis inducing ligand,sTRAIL ELISA KIT	Human	TNFSF10	3.9 pg/ml	15.6 pg/ml - 1000pg/ml
CSB-E07371m	Mouse soluble tumor necrosis factor-related apoptosis inducing ligand,sTRAIL ELISA KIT	Mouse	TNFSF10	7.8 pg/ml	31.25 pg/ml - 2000pg/ml
CSB-E16492h	Human soluble tumor necrosis factor-like weak inducer of apoptosis (sTWEAK)ELISA Kit	Human	TNFSF12	0.078 ng/ml	0.312 ng/ml-20 ng/ml
CSB-EL023987MO	Mouse Tumor necrosis factor ligand superfamily member 12(TNFSF12) ELISA kit	Mouse	TNFSF12	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-EL023991HU	Human Tumor necrosis factor ligand superfamily member 14(TNFSF14) ELISA kit	Human	TNFSF14	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E07380r	Rat Tumor necrosis factor soluble receptor II ,TNFsR- II ELISA KIT	Rat	TNFSR2	0.156 ng/mL	0.625 ng/mL-40 ng/mL
CSB-EQ023955HU	Human Tumor necrosis factor α (TNF-α) antibody ELISA kit	Human	TNF-α Ab	Request Information	Request Information
CSB-EL024016MO	Mouse Troponin T, cardiac muscle(TNNT2) ELISA kit	Mouse	TNNT2	7.8 pg/ml	31.25 pg/ml-2000pg/ml
CSB-E04748h	Human tumor necrosis factor-related apoptosis-inducing ligand receptor 3,TRAIL-R3 ELISA Kit	Human	TRAIL-R3	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-E09623h	Human Trypsin-2(PRSS2) ELISA kit	Human	Try	0.19 ng/ml	0.78 ng/ml-50 ng/ml
CSB-E09316h	Human thymic stromal lymphopoietin,TSLP ELISA Kit	Human	TSLP	0.078 ng/ml	0.312 ng/ml-20 ng/ml
CSB-EL025383MO	Mouse Thioredoxin-interacting protein(TXNIP) ELISA kit	Mouse	TXNIP	0.11 ng/ml	0.45 ng/ml-30 ng/ml
CSB-E11718h	Human Vascular Endothelial cell Growth Factor(VEGF) ELISA KIT	Human	VEGFA	25.297 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E04756m	Mouse Vascular Endothelial cell Growth Factor,VEGF ELISA KIT	Mouse	VEGFA	0.857 pg/mL	3.906 pg/mL-250 pg/mL
CSB-E04757r	Rat Vascular Endothelial cell Growth Factor,VEGF ELISA KIT	Rat	VEGFA	0.97 pg/mL	3.9 pg/mL-250 pg/mL
CSB-E04760h	Human Vascular Endothelial cell Growth Factor D,VEGF-D ELISA KIT	Human	VEGFD	7.81 pg/ml	31.25 pg/ml - 2000 pg/ml
CSB-E04762m	Mouse Vascular endothelial cell growth factor receptor 1,VEGFR-1/Flt1 ELISA kit	Mouse	VEGFR-1/Flt1	0.39 ng/mL	1.56 ng/mL-100 ng/mL
CSB-E04763h	Human Vascular endothelial cell growth factor receptor 2,VEGFR-2/Flk-1 ELISA kit	Human	VEGFR-2/Flk-1	0.098 ng/mL	0.39 ng/mL-25 ng/mL
CSB-E04764m	Mouse Vascular endothelial cell growth factor receptor 2,VEGFR-2/Flk-1 ELISA kit	Mouse	VEGFR-2/Flk-1	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E04765h	Human Vascular endothelial cell growth factor receptor 3,VEGFR-3/Flt-4 ELISA Kit	Human	VEGFR-3/Flt-4	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E04766m	Mouse Vascular endothelial cell growth factor receptor 3 (VEGFR-3/Flt-4) ELISA kit	Mouse	VEGFR-3/Flt-4	19.5 pg/ml	78 pg/ml-5000 pg/ml
CSB-E08712h	Human lymphotactin,Lptn/LTN ELISA Kit	Human	XCL1	0.0975 ng/ml	0.39 ng/ml - 25 ng/ml
CSB-E08714m	Mouse lymphotactin,Lptn/LTN ELISA Kit	Mouse	XCL1	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E04533h	Human eosinophil chemotactic factor,ECF ELISA Kit	Homo sapiens (Human)	CCL11	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E07997h	Human macrophage inflammatory protein 5,MIP-5 ELISA Kit	Homo sapiens (Human)	CCL15	7.81 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E04670h	Human Macrophage Inflammatory Protein-3β,MIP-3β ELISA kit	Homo sapiens (Human)	CCL19	9.75 pg/ml	39 pg/ml-2500 pg/ml
CSB-E04655h	Human monocyte chemotactic protein 1/monocyte chemotactic and activating factor,MCP-1/MCAF ELISA kit	Homo sapiens (Human)	CCL2	1.56 pg/ml	6.25 pg/ml-400 pg/ml
CSB-E07430m	Mouse monocyte chemotactic protein 1/monocyte chemotactic and activating factor,MCP-1/MCAF ELISA kit	Mus musculus (Mouse)	CCL2	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-E04667h	Human Macrophage Inflammatory Protein-3α,MIP-3α ELISA kit	Homo sapiens (Human)	CCL20	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E07309h	Human Myeloid Progenitor Inhibitory Factor 1,MPIF-1 ELISA Kit	Homo sapiens (Human)	CCL23	9.7 pg/mL	39 pg/mL-2500 pg/mL
CSB-E10360h	Human Eotaxin 2/CCL24 ELISA Kit	Homo sapiens (Human)	CCL24	67.660 pg/mL	78.125 pg/mL-5000 pg/mL
CSB-E10019h	Human B-Lymphocyte Chemoattractant 1,BLC-1 ELISA Kit	Homo sapiens (Human)	CXCL13	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-E07419r	Rat macrophage inflammatory protein 2,MIP-2 ELISA kit	Rattus norvegicus (Rat)	CXCL2	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-E08178h	Human Epithelial neutrophil activating peptide 78 (ENA-78/CXCL5) ELISA Kit	Homo sapiens (Human)	CXCL5	0.567 ng/mL	0.625 ng/mL-40 ng/mL
CSB-E08179r	Rat Epithelial neutrophil activating peptide 78,ENA-78 ELISA Kit	Rattus norvegicus (Rat)	CXCL5	0.078 ng/mL	0.312 ng/mL-20 ng/mL
CSB-E08180m	Mouse Epithelial neutrophil activating peptide 78,ENA-78 ELISA Kit	Mus musculus (Mouse)	CXCL5	15.6 pg/ml	62.5 pg/ml - 4000 pg/ml
CSB-E09990h	Human granulocyte chemotactic protein-2,GCP-2 ELISA Kit	Homo sapiens (Human)	CXCL6	19.5 pg/mL	78 pg/mL-5000 pg/mL
CSB-E04539m	Mouse Erythropoietin,EPO ELISA Kit	Mus musculus (Mouse)	EPO	7.8 pg/mL	31.25 pg/mL-2000 pg/mL
CSB-E09027h	Human Heparin-Platelet factor4 Complex antibody(IgG)ELISA Kit	Homo sapiens (Human)	HPF4 Complex Ab (IgG)	Request Information	Request Information
CSB-EQ027171HU	Human anti-IL1 autoantibody ELISA kit	Homo sapiens (Human)	IL1 auto-Ab	Request Information	Request Information
CSB-E13789h	Human IL-1 Receptor Like 1(IL1RL1)ELISA Kit	Homo sapiens (Human)	IL1RL1	0.078 ng/mL	0.312 ng/mL-20 ng/mL
CSB-EQ027170HU	Human anti-IL6 autoantibody ELISA kit	Homo sapiens (Human)	IL6 auto-Ab	Request Information	Request Information
CSB-E04572m	Mouse glucoprotein 130,gp130 ELISA KIT	Mus musculus (Mouse)	IL6ST	19.5 pg/ml	78 pg/ml - 5000 pg/ml
CSB-E08721h	Human myeloperoxidase,MPO ELISA Kit	Homo sapiens (Human)	MPO	0.039 ng/mL	0.156 ng/mL-10 ng/mL
CSB-E07882h	Human Platelet Factor 4,PF-4 ELISA Kit	Homo sapiens (Human)	PF4	1.56 ng/mL	6.25 ng/mL-400 ng/mL
CSB-E07884m	Mouse Platelet Factor 4,PF-4 ELISA Kit	Mus musculus (Mouse)	PF4	1.603 ng/mL	3.125 ng/mL-200 ng/mL
CSB-E04692h	Human Osteoprotegerin,OPG ELISA KIT	Homo sapiens (Human)	TNFRSF11B	1.95 pg/mL	7.8 pg/mL-500 pg/mL
CSB-E07404r	Rat Osteoprotegerin,OPG ELISA KIT	Rattus norvegicus (Rat)	TNFRSF11B	0.195 ng/ml	0.78 ng/ml-50 ng/ml
CSB-E04693m	Mouse Osteoprotegerin,OPG ELISA KIT	Mus musculus (Mouse)	TNFRSF11B	3.9 pg/mL	15.6 pg/mL-1000 pg/mL
CSB-E04525m	Mouse Cluster of differentiation 30,CD30 ELISA Kit	Mus musculus (Mouse)	TNFRSF8	3.9 pg/mL	15.6 pg/mL-1000 pg/mL

CSB-E05127m	Mouse soluble receptor activator of nuclear factor- κ B ligand.sRANKL ELISA Kit	Mus musculus (Mouse)	TNFSF11	3.125 pg/mL	12.5 pg/mL-800 pg/mL
CSB-E04715m	Mouse Soluble Cluster of differentiation 30 ligand.sCD30L ELISA Kit	Mus musculus (Mouse)	TNFSF8	0.039 ng/ml	0.156 ng/ml - 10 ng/ml
CSB-E04577h-IS	Human Interferon γ (IFN- γ) ELISA Kit	Homo sapiens (Human)		1.95 pg/ml	7.8 pg/ml-500 pg/ml
CSB-E04593h-IS	Human Interleukin 10,IL-10 ELISA Kit	Homo sapiens (Human)	IL10	0.975 pg/ml	3.9 pg/ml-250 pg/ml
CSB-E04599h-IS	Human Interleukin 12,IL-12/P70 ELISA Kit	Homo sapiens (Human)	IL-12/P70	0.975 pg/ml	3.9 pg/ml-250 pg/ml
CSB-E04608m-IS	Mouse Interleukin 17,IL-17 ELISA Kit	Mus musculus (Mouse)	IL17A	2.735 pg/ml	10.938 pg/ml-700 pg/ml
CSB-E04627m-IS	Mouse Interleukin 2,IL-2 ELISA kit	Mus musculus (Mouse)	IL2	0.975 pg/ml	3.9 pg/ml-250 pg/ml
CSB-E04639m-IS	Mouse Interleukin 6,IL-6 ELISA Kit	Mus musculus (Mouse)	IL6	0.975 pg/ml	3.9 pg/ml-250 pg/ml
CSB-E04741m-IS	Mouse Tumor necrosis factor α ,TNF- α ELISA Kit	Mus musculus (Mouse)	TNF	1.95 pg/ml	7.8 pg/ml-500 pg/ml
CSB-E08053h-IS	Human Interleukin 1 β ,IL-1 β ELISA Kit	Homo sapiens (Human)	IL1B	0.586 pg/ml	2.344 pg/ml-150 pg/ml
CSB-E08054m-IS	Mouse Interleukin 1 β ,IL-1 β ELISA Kit	Mus musculus (Mouse)	IL1B	3.125 pg/ml	12.5 pg/ml-800 pg/ml
CSB-E04725h-IS	Human Transforming Growth factor β 1,TGF- β 1 ELISA kit	Homo sapiens (Human)	TGFB1	5.8 pg/ml	23.5 pg/ml-1500 pg/ml