

Competent cells selection guide

	Blue-white capable	Phage T1 resistant	Restriction-deficient (End A-/RecA)	Clone methylated (RMS)	Inducible plasmid copy number	Fepisome	Antibiotic resistance	Electrocompetent: efficiency	Chemically competent: efficiency	
Routine cloning and library construction										
E. cloni® 10G (Classic)	✓	✓	✓	✓	✓	Ø	Ø	Str	≥5 × 10 ⁹	10 ⁶ –10 ⁹
E. cloni 10G (Elite)	✓	✓	✓	✓	✓	Ø	Ø	Str	≥2 × 10 ¹⁰	–
E. cloni 10G (Supreme)	✓	✓	✓	✓	✓	Ø	Ø	Str	≥4 × 10 ¹⁰	–
E. cloni 10G F' (Elite)	✓	✓	✓	✓	✓	Ø	✓	Str	≥2 × 10 ¹⁰	–
E. cloni 5-alpha	✓	✓	✓	Ø	Ø	Ø	Ø	–	–	≥1 × 10 ⁸
Unstable inserts or toxic gene products										
Copycutter™ EPI400™	✓	✓	✓	✓	✓	✓	Ø	Str	≥1 × 10 ¹⁰	≥1 × 10 ⁷
Large inserts, plasmids, or fosmids										
TransforMax™ EC100™	✓	Ø	✓	✓	✓	Ø	Ø	Str	≥1 × 10 ¹⁰	≥5 × 10 ⁸
TransforMax EPI300™	✓	Ø	✓	✓	✓	a	Ø	Str	≥1 × 10 ¹⁰	≥5 × 10 ⁸
TransforMax EPI300™-T1R	✓	✓	✓	✓	✓	a	Ø	Str	≥1 × 10 ¹⁰	–
BAC cloning										
BAC-Optimised Replicator™ v2.0	✓	Ø	✓	✓	✓	a	Ø	Amp, Str	≥1 × 10 ¹⁰	–
E. cloni BAC-Optimised 10G	✓	Ø	✓	✓	✓	Ø	Ø	Str	≥1 × 10 ¹⁰	–
BigEasy™-TSA	b	✓	✓	✓	✓	c	Ø	Amp, Str	≥2 × 10 ¹⁰	–
Phage display										
TG1	✓	Ø	Ø	Ø	d	Ø	✓	–	≥4 × 10 ¹⁰	–
ER2738	✓	✓	Ø	Ø	d	Ø	✓	Tet	≥2 × 10 ¹⁰	–
SS320 (MC1061 F')	✓	Ø	Ø	Ø	e	Ø	✓	Tet, Str	≥4 × 10 ¹⁰	–
MC1061 F-	Ø	Ø	Ø	Ø	e	Ø	Ø	Str	≥4 × 10 ¹⁰	–
CRISPR sgRNA libraries and lentiviral plasmid cloning										
Endura™	Ø	Ø	f	Ø	g	Ø	Ø	Str	≥1 × 10 ¹⁰	≥1 × 10 ⁷
Site-directed mutagenesis										
CJ236	Ø	Ø	Ø	Ø	Ø	Ø	✓	Cam	≥1 × 10 ⁹	–

Amp, ampicillin; Cam, chloramphenicol; Str, streptomycin; Tet, tetracycline

a: Requires OriV vectors (e.g., CopyRight® and pCCFOS)

b: Requires IPTG

c: Only with pJAZZ® vectors

d: mcrB– only; not suitable for cloning eukaryotic genomic DNA

e: mcrA– mcrB–; not suitable for cloning eukaryotic genomic DNA

f: RecA13 only

g: mcrB– mrr–; not suitable for cloning eukaryotic genomic DNA

	Primary promoter	Background expression control	Low endotoxin	Antibiotic resistance	Electrocompetent: efficiency	Chemically competent: efficiency
Protein expression: routine						
E. cloni EXPRESS BL21(DE3)	T7	Ø	Ø	–	$\geq 5 \times 10^9$	$\geq 1 \times 10^7$
HI-Control™ BL21(DE3)	T7	✓	Ø	Gen	–	$\geq 1 \times 10^7$
HI-Control™ 10G	P_{lac}^+ P_{lac}^+ P_{trc}^+ $T5_{lac}$ $rhaP_{BAD}$	✓	Ø	Gen, Str	–	$\geq 1 \times 10^9$
ClearColi® BL21(DE3)	T7	Ø	✓	–	$\geq 1 \times 10^9$	–
Protein expression: toxic products						
OverExpress™ C41(DE3)	T7	Ø	Ø	–	$\geq 1 \times 10^9$	$\geq 1 \times 10^6$
OverExpress C41(DE3) pLysS	T7	✓	Ø	Cam	–	$\geq 1 \times 10^6$
OverExpress C43(DE3)	T7	Ø	Ø	–	$\geq 1 \times 10^9$	$\geq 1 \times 10^6$
OverExpress C43(DE3) pLysS	T7	✓	Ø	Cam	–	$\geq 1 \times 10^6$

Cam, chloramphenicol; Gen, gentamicin; Str, streptomycin; rhaPBAD, rhamnose-inducible promoter

Not finding a strain suitable for your application? Do you need a different size or bulk purchase? Lucigen's Custom Competent Cell Preparation Services provide:

- **Support** - Personalised service from initial discussion through delivery
- **Flexibility** - Choose from the following options to match your specific needs:
 - Style:** Chemically competent or electrocompetent cells
 - Dispense volume:** bulk or small aliquots
 - Format:** 96-well plates or tubes of your choice
 - Strains:** Lucigen or any BSL1 E. coli strain of choice that meets your application needs
- **Transformation efficiency** - High efficiencies that match your needs, up to an industry-leading 4×10^{10} cfu/µg DNA
- **Rapid turnaround time** - Average 2–3 weeks
- **Quality** - ISO 13485 Certified with 10+ years of manufacturing excellence, available for OEM and partnerships

Request a free, no-obligation quote today!

Visit lucigen.com/customcells or e-mail: bizdev@lucigen.com

Integrated tools. Accelerated science.




 @LGCBiosearch | lucigen.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2018. All rights reserved. GEN/0526/SW/1118 DS-181129.01

For research or investigational use only.



**LUBIO
SCIENCE**

Your distributor in Switzerland

LubioScience GmbH
Baumackerstrasse 24
8050 Zürich
+41 (0)41 417 02 80

info@lubio.ch
www.lubio.ch

**BIOSEARCH™
TECHNOLOGIES**

GENOMIC ANALYSIS BY LGC